



InertSearch™ for GC

InertCap® Applications

Organic Solvents

Data No. GA202-0410



InertSearch™ for GC

InertCap® Applications

Organic Solvents

Data No. GA202-0410

Conditions

System	: GC/FID
Column	: TC-BOND Q 0.53 mm I.D. x 30 m df = 20 µm
Col. Cat. No.	: 1010-40243
Col. Temp.	: 40 °C - 5 °C/min - 250 °C
Carrier Gas	: He 30 kPa
Injection	: Split flow 150 mL/min 250 °C
Detection	: FID Range 10 ¹ 250 °C
Sample Size	: Mixed evenly 0.2 µL

Analyte

1. Methanol
2. Ethanol
3. Acetonitrile
4. Acetone
5. Dichloromethane
6. 2-Propanol(Isopropyl alcohol)
7. Methyl acetate
8. Ethyl ether
9. 1-Propanol
10. *tert*-Butanol
11. Tetrahydrofuran
12. Methyl ethyl ketone(MEK)
13. Chloroform
14. Ethyl acetate
15. 2-Butanol
16. 1,2-Dichloroethane
17. Methyl cellosolve
18. 2-Methyl-1-propanol(Isobutyl alcohol)
19. *n*-Hexane
20. 1,1,1-Trichloroethane
21. Carbon tetrachloride
22. Benzene
23. 1-Butanol
24. Trichloroethylene
25. 1,4-Dioxane
26. Isopropyl acetate