


FAMES (PUFA, animal source)
MXT®-WAX

30m, 0.28mm ID, 0.25µm MXT®-WAX (cat.# 70624)
0.1µL split injection of PUFA 2 mix.

Oven temp.: 160°C to 250°C @ 2°C/min. (hold 10 min.)
Inj. & det. temp.: 260°C
Carrier gas: hydrogen
Linear velocity: 40cm/sec.
FID sensitivity: 8 x 10⁻¹¹ AFS
Split ratio: 20:1



1. C14:0
2. C16:0
3. C16:1n7
4. C18:0
5. C18:1n9
6. C18:1n7
7. C18:2n6
8. C18:3n3
9. C18:4n3
10. C20:1n9
11. C20:2n6
12. C20:3n6
13. C20:4n6
14. C20:5n3
15. C22:4n6
16. C22:5n3
17. C22:6n3
18. C24:1n9