

# Alcoholic Beverage Standard: Acids and Esters

## Stabilwax®-DA

Stabilwax®-DA 30m, 0.18mm ID, 0.18 $\mu$ m (cat.# 550752)  
Inj.: 1 $\mu$ L splitless (hold 0.5 min.) at conc. shown in  
peak list, in ethyl acetate, 4mm ID splitless liner  
w/wool (cat.# 20814-202.1)  
Inj. temp.: 240°C  
Carrier gas: hydrogen  
Make-up gas: nitrogen  
Linear velocity: 28psi @ 240°C  
Oven temp.: 70°C to 240°C at 12°C/min. (hold 3 min.)  
Det.: FID

