


## Triglyceride distribution in olive oil

**BPX5**

Column part number	054133		
Phase	BPX5	Final temperature	360°C, 5 min
Column	12 m x 0.53 x 0.25 µm	Detector	FID
Initial temperature	100°C, 0.5 min	Carrier gas	He, 6 psi
Rate 1	50°C/min to 280°C	Injection mode	On-column
Rate 2	10°C/min to 360°C	Injection volume	1 µL



For more information visit [www.trajanscimed.com](http://www.trajanscimed.com) or contact [techsupport@trajanscimed.com](mailto:techsupport@trajanscimed.com)